


Doing A Bypass Procedure?

ViSiGi® is integral for your gastric bypass, aiding in four key steps with just one insertion for increased efficiency.

- Decompression
- Pouch Creation
- Sizing of the Gastrojejunostomy
- Leak Testing



Short Hole
Pattern



ViSiGi®
for Bypass
Procedures

Standard
Hole Pattern



ViSiGi®
for Sleeve
Gastrectomy

ViSiGi® for Bypass Features

- **Selection Valve**

Vented close position allows for safe adjustment of the device. Open position allows for decompression and internal stabilization.



- **Safe Suction**

Integral suction control automatically reduces High-Vacuum to clinically effective levels to ensure maximum delineation.



- **Short Hole Pattern**

A shorter hole pattern allows you to size and leak-test your gastrojejunostomy.

P/N	Description	Ordering
5232S	ViSiGi 3D® for Bypass, 32Fr	1 Box (5 Systems)
5236S	ViSiGi 3D® for Bypass, 36Fr	1 Box (5 Systems)
5240S	ViSiGi 3D® for Bypass, 40Fr	1 Box (5 Systems)
5232SB	ViSiGi 3D® for Bypass with Bulb, 32Fr	1 Box (5 Systems)
5236SB	ViSiGi 3D® for Bypass with Bulb, 36Fr	1 Box (5 Systems)
5240SB	ViSiGi 3D® for Bypass with Bulb, 40Fr	1 Box (5 Systems)